

# Quick Cut

Model: C0046

*a measurable difference...*

**IDM**<sup>®</sup>  
instruments

Quick Cut is an easy to operate apparatus used for cutting of corrugated cardboard for sample preparation.

## Applications:

- Corrugated Cardboard

## Features:

- Easy operation
- Engraved Cutting Guides
- 5 spare blades
- Thickness: Up to 10mm

## Benefits:

- Easy to use
- Fast results
- Accurate

## Standards:

- N/A

## Options:

- Different sized Cutter available upon request

## Connections:

- N/A

## Dimensions:

- **H:** 80mm
- **W:** 400mm
- **D:** 480mm
- **Weight:** 10kg

