

# Compression Tester-65 Tonne

Model: C0011

*a measurable difference...*

**IDM**<sup>®</sup>  
instruments

The 65 Tonne Press is easy to use with an electric hydraulic pump with a vertical above mounted cylinder. Designed for various compression tests it comes fitted with load cells to measure the load applied. Having a large work area, this machine is adaptable to many applications.

## Applications:

- Various

## Features:

- Hydraulic Power Unit
- Electronic Load Cells with an accuracy of +/- 1% FS
- Position Measuring Device

## Benefits:

- Easy to use
- Fast results
- Accurate

## Standards:

- N/A

## Options:

- N/A

## Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or  
110 VAC @ 60 HZ  
(please specify when ordering)

## Dimensions:

- **H:** TBA
- **W:** TBA
- **D:** TBA
- **Weight:** TBA

